

# Ultrafast Dental Laser

## Wireless is Effortless

Amazing Laser for Everyday Practice

*No More Setup Time and Messy Cords*

3W Continuous Power

5W Pulsed Power



**Make it Ultrafast**

**Hemostasis**

**Frenectomy**

**Gingivectomy**

**Incision & Excision**

**Root Canal Therapy**

**Canker Sores / Aphthous Ulcers**

**Periodontal Sulcular Debridement**

**Gingival Troughing for Impression**

# Ultrafast Dental Laser

## Amazing Power at Your Fingertips

### Wireless Powerful Design

- Up to 3W CW Power for ultrafast soft tissue sculpting
- Up to 5W Pulsed Power for precision cutting
- 810nm wavelength for advanced performance

### Making it Easy

- Distinctive green aiming beam - higher contrast making it easier to see the treatment spot
- Unique Laser Shield
  - Automatic Eye Protection for Operators (including Loupe)
- Disposable Fiber tips are:
  - o Time-saving
  - o Self-bendable
  - o Avoid cross-contamination

### Super Slim and Ergonomic

- Only weighs 3.5oz
- Measures at 6.4 inch

### Built to Last

Aircraft aluminum durable construction distributes heat, prevent scratches, discoloration and damage from accidental falls.

### Smart Lithium Battery

- 45+ minutes run time at nominal power
- 90 minutes fast recharge time
- Constant output throughout battery run time
- Charges directly or through charging stand
- Compatible with FUSION curing lights

